

Nifuratel

Catalog No: #S1526


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nifuratel
Brief Description	Inhibitors
Purification	99.00%
Target Name	Antibacterial inhibitor
Calculated MW	285.28
Formulation	C10H11N3O5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CSCC1CN(C(=O)O1)/N=C/C2=CC=C(O2)[N+](=O)[O-] Pathways: Microbiology&virology Receptor: Antibacterial Boling
pt: Melting pt: 176-178°C Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Polatti F. Curr Clin Pharmacol. 2012 Feb 1;7(1):36-40.

Note: This product is for in vitro research use only