

Pidotimod

Catalog No: #S1532


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pidotimod
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	244.27
Formulation	C ₉ H ₁₂ N ₂ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: C1CC(=O)N[C@@H]1C(=O)N2CSC[C@H]2C(=O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: 663C
 Melting pt: 194-198C
 Solubility: Water: 10 mg/mL (40.93 mM); DMSO: 49 mg/mL (200.59 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Clemente E, et al. *Arzneimittelforschung*, 1994, 44(12A), 1490-1494.

Note: This product is for in vitro research use only