

## Valganciclovir hydrochloride

Catalog No: #S1533


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Valganciclovir hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	390.82
Formulation	C <sub>14</sub> H <sub>23</sub> CIN <sub>6</sub> O <sub>5</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: CC(C)[C@@H](C(=O)OCC(CO)OCN1C=NC2=C1NC(=NC2=O)N)N  
 CIP Pathways: Others  
 Receptor: CMV  
 Boiling pt: Melting pt: 162-164°C  
 Solubility: Water: 78 mg/mL (199.58 mM); DMSO: 78 mg/mL (199.58 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Sugawara M, et al. J Pharm Sci, 2000, 89(6), 781-789.

Note: This product is for in vitro research use only