

Darifenacin HBr

Catalog No: #S1534


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Darifenacin HBr
Brief Description	Inhibitors
Purification	99.50%
Target Name	AChR antagonist
Calculated MW	505.45
Formulation	C28H31BrN2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:C1CN(C[C@@H]1C(C2=CC=CC=C2)(C3=CC=CC=C3)C(=O)N)CCC4=CC5=C(C=C4)OCC5.BrPathways:NeuroscienceReceptor:M3
 mACHRBoling pt:614.3CMelting pt:228-230CSolubility:DMSO: 117 mg/mL (230.56 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Moriya H, et al. Life Sci. 1999; 64(25):2351-8.

Note: This product is for in vitro research use only