


Description

Product Name	Promestriene
Brief Description	Inhibitors
Purification	98.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	328.49
Formulation	C22H32O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CCCOC1=CC2=C(C=C1)[C@H]3CC[C@]4([C@H]([C@@H]3CC2)CC[C@@H]4OC)CPathways:Endocrinology/HormonesReceptor:ERBoiling pt:436.8CMelting pt:64-66CSolubility:Water: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Del Pup L, et al. AntiY Drugs. 2013 Nov;24(10):989-98.

Note: This product is for in vitro research use only