

## L-Cysteine methyl ester hydrochloride(Acdrile)

Catalog No: #S1547


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |  |
|-------------------|--|
| Product Name      | L-Cysteine methyl ester hydrochloride(Acdrile) |
| Brief Description | Inhibitors                                     |
| Purification      | 98.00%   |
| Target Name       | Amino Acids and Derivatives inhibitor          |
| Calculated MW     | 171.65   |
| Formulation       | C4H10ClNO2S                                    |
| Storage           | 3 years -20°C powder;2 years -80°C in solvent; |

## Images



## Product Description

Research Area:OtherSMILES:COC(=O)[C@H](CS)N.ClPathways:OthersReceptor:Amino Acids and DerivativesBoiling pt:197.2CMelting pt:142CSolubility:Water: 10 mMAppearance:White Needle CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. CHARPENTIER, DARROUZET. Rev Laryngol Otol Rhinol (Bord). 1956 Jan-Feb;77(1-2):147-53.

Note: This product is for in vitro research use only