

Tegaserod maleate

Catalog No: #S1551


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tegaserod maleate
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist
Calculated MW	417.46
Formulation	C20H27N5O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CCCCCN=C(N)NN/C=C/1C=NC2=C1C=C(C=C2)OC.C(=C\C(=O)O)\C(=O)OPathways:GPCR/G Protein;
NeuroscienceReceptor:5-HT; Sert (Sodium-dependent)Boiling pt:Melting pt:180-183CSolubility:DMSO: 35mg/mLAppearance:White
PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be
stored below -20 °C for several month.

References

1. Beattie DT, et al. Br J Pharmacol. 2004 Nov;143(5):549-60.

Note: This product is for in vitro research use only