

Iohexol


Catalog No: #S1554

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Iohexol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	821.14
Formulation	C19H26I3N3O9
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC(=O)N(CC(O)CO)C1=C(I)C(C(=O)NCC(O)CO)=C(I)C(C(=O)NCC(O)CO)=C1I Pathways: Others Receptor: Others Boiling pt: 891.5°C Melting pt: 254-256°C Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yin Q, et al. J Am Chem Soc. 2013 Sep 18; 135(37):13620-3.

Note: This product is for in vitro research use only