

Ornidazole

Catalog No: #S1556


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ornidazole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	219.63
Formulation	C7H10ClN3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CC1=NC=C(N1CC(O)CCl)[N+](=O)[O-]
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 85-90°C
Solubility: Ethanol: 44 mg/mL (200.33 mM); Water: 2 mg/mL (9.1 mM); DMSO: 44 mg/mL (200.33 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rutgeerts P et al. Gastroenterology, 2005, 128(4):856-861.

Note: This product is for in vitro research use only