

Rebamipide


Catalog No: #S1562

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rebamipide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Free radical scavengers inhibitor
Calculated MW	370.79
Formulation	C19H15ClN2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: OC(=O)C(CC1=CC(=O)NC2=CC=CC=C12)NC(=O)C1=CC=C(Cl)C=C1 Pathways: oxidation-reduction Receptor: Free radical scavengers
Boiling pt: 695°C Melting pt: 288-290°C Solubility: DMSO: 74 mg/mL (199.57 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Matsuda T, et al. Drugs R D. 2003;4(1):19-28.

Note: This product is for in vitro research use only