

EphA7 Polyclonal Antibody

Catalog No: #40888

Package Size: #40888-1 50ul #40888-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EphA7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	EphA7 Polyclonal Antibody detects endogenous levels of EphA7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human EphA7.
Target Name	EphA7
Other Names	EPHA7; EHK3; HEK11; Ephrin type-A receptor 7; EPH homology kinase 3; EHK-3; EPH-like kinase 11; EK11; hEK11
Accession No.	Swiss-Prot: Q15375NCBI Gene ID: 2045
Uniprot	Q15375
GeneID	2045;
SDS-PAGE MW	112kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T-UV cells using EphA7 Polyclonal Antibody

Note: This product is for in vitro research use only