

Cisplatin

Catalog No: #S1564

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cisplatin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Caspase activator
Calculated MW	300.05
Formulation	H6Cl2N2Pt
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:Cl[Pt+ +]ClPathways:Apoptosis; Proteases/ProteasomeReceptor:Caspase-3; Caspase-9Boiling pt:157CMelting pt:77-78CSolubility:DMF: 12 mg/mL (39.99 mM); 12 mg/mL (39.99 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Siddik ZH. Oncogene, 2003, 22(47), 7265-7279.

Note: This product is for in vitro research use only