

Flunisolide

Catalog No: #S1573


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flunisolide
Brief Description	Inhibitors
Purification	98.00%
Target Name	GR inhibitor
Calculated MW	434.5
Formulation	C ₂₄ H ₃₁ FO ₆
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12C[C@@H]([C@H]3[C@H]([C@@H]1C[C@@H]4[C@]2(OC(O4)(C)C(=O)CO)C[C@@H](C5=CC(=O)C=C[C@]35C)F)O

Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Hochhaus G, Moellmann HW. Agents Actions. 1990 Jun;30(3-4):377-80.

Note: This product is for in vitro research use only