

## lopamidol

Catalog No: #S1576


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	lopamidol
Brief Description	Inhibitors
Purification	99.47%
Target Name	Others
Calculated MW	777.09
Formulation	C17H22I3N3O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C[C@@H](C(=O)NC1=C(C(=C(C(=C1)C(=O)NC(CO)CO))C(=O)NC(CO)CO))O Pathways: Others Receptor: Others Boling  
 pt: Melting pt: >320C Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Weiland FL, et al. Acta Radiol. 2014 Jul;55(6):715-24.

Note: This product is for in vitro research use only