

Midodrine hydrochloride

Catalog No: #S1577


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Midodrine hydrochloride
Brief Description	Inhibitors
Purification	99.36%
Target Name	Adrenergic Receptor agonist
Calculated MW	290.74
Formulation	C ₁₂ H ₁₉ CIN ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: Cl.COC1=CC(C(O)CNC(=O)CN)=C(OC)C=C1
 Pathways: GPCR/G Protein Receptor: α-adrenergic receptor
 Boiling pt: Melting pt: 192-193°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Buckner SA, et al. Eur J Pharmacol. 2002 Aug 2;449(1-2):159-65.

Note: This product is for in vitro research use only