

Pargyline hydrochloride

Catalog No: #S1578


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pargyline hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	MAO inhibitor
Calculated MW	195.69
Formulation	C ₁₁ H ₁₄ NCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: Cl.CN(CC#C)CC1=CC=CC=C1
Pathways: Metabolism
Receptor: MAO
Boiling pt: Melting pt: 160-163°C
Solubility: DMSO: 38 mg/mL
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Patkar AA, et al. CNS Spectr. 2006, 11(5):363-75.

Note: This product is for in vitro research use only