

Mestranol

Catalog No: #S1579


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mestranol
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	310.43
Formulation	C ₂₁ H ₂₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@]2(C#C)O)CCC4=C3C=CC(=C4)OCPathways:Endocrinology/HormonesReceptor:ERBo
 ling pt:Melting pt:153-155CSolubility:Ethanol: 25 mg/mL (80.53 mM); DMSO: 15 mg/mL (48.32 mM)Appearance:White PowderRemark:For obtaining a
 higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several
 month.

References

1. IARC Monogr Eval Carcinog Risk Chem Hum. 1979 Dec;21:257-78.

Note: This product is for in vitro research use only