

## Naproxen

Catalog No: #S1582


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Naproxen
Brief Description	Inhibitors
Purification	99.60%
Target Name	COX inhibitor
Calculated MW	230.27
Formulation	C14H14O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:C[C@@H](C1=CC2=C(C=C1)C=C(C=C2)OC)C(=O)O

Pathways:Immunology/Inflammation;

NeuroscienceReceptor:COXBoling pt:Melting pt:152-154CSolubility:Soluble in mEthanolAppearance:White PowderRemark:For obtaining a higher

solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Albertini R, et al. Inflamm Res. 2007 Jun;56(6):228-9.

Note: This product is for in vitro research use only