

Tenoxicam

Catalog No: #S1586


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tenoxicam
Brief Description	Inhibitors
Purification	98.00%
Target Name	COX inhibitor
Calculated MW	337.37
Formulation	C13H11N3O4S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CN1C(C(=O)NC2=CC=CC=N2)=C(O)C2=C(C=CS2)S1(=O)=OPathways:Immunology/Inflammation;
NeuroscienceReceptor:COXBoling pt:Melting pt:209-213CSolubility:DMSO: 68 mg/mL (201.55 mM)Appearance:Yellow PowderRemark:For obtaining
a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several
month.

References

1. Yamada M, et al. J Pharmacol Exp Ther. 1997 May;281(2):1005-12.

Note: This product is for in vitro research use only