

Exemestane

Catalog No: #S1587


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Exemestane
Brief Description	Inhibitors
Purification	99.35%
Target Name	Aromatase inhibitor
Calculated MW	296.4
Formulation	C ₂₀ H ₂₄ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: C[C@]12CC[C@H]3[C@H]([C@@H]1CCC2=O)CC(=C)C4=CC(=O)C=C[C@]34 Pathways: Endocrinology/Hormones Receptor: Aromatase (human); Aromatase (rat) Boiling pt: Melting pt: 155.13°C Solubility: Ethanol: 15 mg/mL (50.6 mM); DMSO: 54 mg/mL (182.18 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wang X, et al. Y Res, 2006, 66(21), 10281-10286.

Note: This product is for in vitro research use only