

## Kemantane( Idramantone)

Catalog No: #S1588


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Kemantane( Idramantone)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others agonist
Calculated MW	166.22
Formulation	C10H14O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC12CC3CC(C1)C(=O)C(C3)C2Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mirzolian NR, et al. Eksp Klin Farmakol. 2013; 76(6): 7-10.

Note: This product is for in vitro research use only