

Nevirapine

Catalog No: #S1595


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nevirapine
Brief Description	Inhibitors
Purification	99.59%
Target Name	HIV Protease inhibitor
Calculated MW	266.3
Formulation	C15H14N4O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CC1=C2NC(=O)C3=C(N=CC=C3)N(C3CC3)C2=NC=C1
 Pathways: Microbiology & virology;
 Proteases/Proteasome Receptor: HIV-1 reverse transcriptase
 Boiling pt: 415.4°C
 Melting pt: 247°C
 Solubility: DMSO: 53 mg/mL (199.02 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Ambrose Z, et al. J Virol. 2009 Apr; 83(8): 3826-33.

Note: This product is for in vitro research use only