

Proadifen hydrochloride

Catalog No: #S1597


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proadifen hydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	P450 inhibitor
Calculated MW	389.96
Formulation	C23H32ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.CCCC(C(=O)OCCN(CC)CC)(C1=CC=CC=C1)C1=CC=CC=C1
 Pathways: Metabolism
 Receptor: CYP
 Boiling pt: 460.8°C
 Melting pt: 122-123°C
 Solubility: Soluble in DMSO
 Appearance: Off-white Color Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Viana AP, Osswald W. *Arzneimittelforschung*. 1970 Jun;20(6):851-3.

Note: This product is for in vitro research use only