

Rimantadine hydrochloride

Catalog No: #S1598


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rimantadine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	215.76
Formulation	C ₁₂ H ₂₁ NHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: Cl.CC(N)C12CC3CC(CC(C3)C1)C2
 Pathways: Others
 Receptor: RNA
 Boling pt: Melting pt: >300C
 Solubility: DMSO: 33 mg/mL; Water: 33 mg/mL; Ethanol 33 mg/mL
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Liu BM, et al. Iran J Pharm Res. 2014 Fall;13(4):1183-90.

Note: This product is for in vitro research use only