

## Ibandronate sodium

Catalog No: #S1599L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ibandronate sodium
Brief Description	Inhibitors
Purification	90.00%
Target Name	HBV inhibitor; Others
Calculated MW	341.21
Formulation	C <sub>9</sub> H <sub>22</sub> NNaO <sub>7</sub> P <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: O.[Na+].[Na+].CCCCN(C)CCC(O)(P(O)(O)=O)P(O)([O-])=O  
 Pathways: Microbiology&virology;  
 Others Receptor: HAP; FPS  
 Boling pt: Melting pt: Solubility: Free soluble in Water. Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Dunford JE, et al. J Pharmacol Exp Ther. 2001 Feb; 296(2):235-42.

Note: This product is for in vitro research use only