

Tiratricol

Catalog No: #S1604


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tiratricol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Thyroid hormone receptor(THR)
Calculated MW	621.93
Formulation	C14H9I3O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:OC(=O)CC1=CC(I)=C(OC2=CC(I)=C(O)C=C2)C(I)=C1Pathways:Endocrinology/HormonesReceptor:Thyroid hormoneBoiling pt:Melting pt:Solubility:Ethanol: 100 mg/mL (160.78 mM); DMSO: 100 mg/mL (160.78 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Carvalho GA, et al. Arq Bras Endocrinol Metabol, 2004, 48(1), 83-92.

Note: This product is for in vitro research use only