

Exherin (trifluoroacetate)


Catalog No: #S1608

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Exherin (trifluoroacetate)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	684.71
Formulation	C ₂₄ H ₃₅ F ₃ N ₃ O ₃ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: C[C@H]1C(=O)N[C@H](C(=O)N[C@@H](CSSC[C@@H](C(=O)N[C@H](C(=O)N1)CC2=CN=CN2)NC(=O)C)C(=O)N)C(C)C.C(=O)C(F)(F)F
 Pathways: Metabolism
 Receptor: N-cadherin
 Boiling pt: Melting pt: Solubility: DMSO: 43 mg/ml
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Perotti A, et al. Ann Oncol. 2009 Apr;20(4):741-5.

Note: This product is for in vitro research use only