

Metaxalone


Catalog No: #S1610

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metaxalone
Brief Description	Inhibitors
Purification	98.84%
Target Name	Others
Calculated MW	221.25
Formulation	C ₁₂ H ₁₅ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CC1=CC(OCC2CNC(=O)O2)=CC(C)=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 121-123°C
Solubility: Water: 91 mg/L . Appearance: White Powder
Remark: For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nirogi RV, et al. J Anal Toxicol. 2006 May;30(4):245-51.

Note: This product is for in vitro research use only