

Hodostin


Catalog No: #S1612

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hodostin
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor
Calculated MW	334.39
Formulation	C13H22N2O6S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: COS(O)(=O)=O
Pathways: Neuroscience
Receptor: AChR
Boiling pt: Melting pt: 175-177°C
Solubility: Water: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Huang TH, et al. Yao Xue Xue Bao. 1982 Jul;17(7):534-9.

Note: This product is for in vitro research use only