

Rosiglitazone maleate

Catalog No: #S1622


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rosiglitazone maleate
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR agonist
Calculated MW	473.12
Formulation	C18H19N3O3SC4H4O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: CN(CCOC1=CC=C(C=C1)CC2C(=O)NC(=O)S2)C3=CC=CC=N3.C(=C\C(=O)O)\C(=O)O
 Pathways: Metabolism Receptor: PPAR γ
 Boiling pt: Melting pt: 235-240°C Solubility: DMSO: 94 mg/mL (198.52 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pftzner A, et al. Vasc Health Risk Manag. 2007; 3(2):211-20.

Note: This product is for in vitro research use only