

## Rivastigmine tartrate

Catalog No: #S1628


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Rivastigmine tartrate
Brief Description	Inhibitors
Purification	99.66%
Target Name	AChR inhibitor
Calculated MW	400.43
Formulation	C <sub>14</sub> H <sub>22</sub> N <sub>2</sub> O <sub>2</sub> C <sub>4</sub> H <sub>6</sub> O <sub>6</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous

systemSMILES:CCN(C)C(=O)OC1=CC=CC(=C1)[C@H](C)N(C)C.[C@@H]([C@H](C(=O)O)O)(C(=O)O)OPathways:NeuroscienceReceptor:AChRBol  
ing pt:316.2CMelting pt:123-125CSolubility:Ethanol: 80 mg/mL (199.79 mM); Water: 80 mg/mL (199.79 mM); DMSO: 42 mg/mL (104.88  
mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a  
while.Stock solution can be stored below -20 °C for several month.

## References

1. Kennedy JS, et al. J Clin Psychopharmacol. 1999 Dec;19(6):513-21.

Note: This product is for in vitro research use only