

Tiopronin

Catalog No: #S1630


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Tiopronin |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 163.19 |
| Formulation | C5H9NO3S |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:ImmunologySMILES:CC(S)C(=O)NCC(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:93-98CSolubility:Ethanol: 32 mg/mL (196.09 mM); Water: 32 mg/mL (196.09 mM); DMSO: 32 mg/mL (196.09 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zheng Z, et al. Ren Fail. 2014 Oct;36(9):1455-60.

Note: This product is for in vitro research use only