

Amitraz

Catalog No: #S1632


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amitraz
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	293.41
Formulation	C19H23N3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system SMILES: CN(C=NC1=C(C)C=C(C)C=C1)C=NC1=C(C)C=C(C)C=C1 Pathways: GPCR/G

Protein Receptor: α -adrenergic receptor Boiling pt: Melting pt: 86-87°C Solubility: DMSO: > 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Del Pino J, et al. Chem Res Toxicol. 2015 Jun 15;28(6):1073-94.

Note: This product is for in vitro research use only