

Epsin 1 Polyclonal Antibody

Catalog No: #40895

Package Size: #40895-1 50ul #40895-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epsin 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Epsin 1 Polyclonal Antibody detects endogenous levels of Epsin 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Epsin 1.
Target Name	Epsin 1
Other Names	EPN1; Epsin-1; EH domain-binding mitotic phosphoprotein; EPS-15-interacting protein 1
Accession No.	Swiss-Prot: Q9Y6I3NCBI Gene ID: 29924
Uniprot	Q9Y6I3
GeneID	29924;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MOuse-brain cells using Epsin 1 Polyclonal Antibody



Western Blot analysis of K562 cells using Epsin 1 Polyclonal Antibody

Note: This product is for in vitro research use only