

Glyburide

Catalog No: #S1634


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glyburide
Brief Description	Inhibitors
Purification	99.77%
Target Name	Potassium Channel inhibitor
Calculated MW	494.01
Formulation	C23H28ClN3O5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: COC1=CC=C(Cl)C=C1C(=O)NCCC1=CC=C(C=C1)S(=O)(=O)NC(=O)NC1CCCCC1 Pathways: Membrane transporter/Ion channel Receptor: Bile salt export pump; Potassium Channel Boiling pt: Melting pt: 168-173°C Solubility: DMSO: 99mg/ml(200.4 mM); Water:

References

1. Pondugula SR, et al. *Physiol Genomics*. 2006 Jan 12;24(2):114-23.

Note: This product is for in vitro research use only