

## N-Acetylneuraminic acid

Catalog No: #S1638


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	N-Acetylneuraminic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	309.27
Formulation	C11H19NO9
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC(=O)N[C@@H]1[C@H](C[C@](O[C@H]1[C@@H]([C@@H](CO)O)O)(C(=O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:184-186CSolubility:Water: 10 mMAppearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Severi E, et al. Microbiology. 2007 Sep;153(Pt 9):2817-22.

Note: This product is for in vitro research use only