

Amlexanox

Catalog No: #S1639


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amlexanox
Brief Description	Inhibitors
Purification	98.00%
Target Name	FGFR inhibitor; IL Receptor antagonist; Others inhibitor
Calculated MW	298.29
Formulation	C16H14N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CC(C)C1=CC2=C(OC3=NC(N)=C(C=C3C2=O)C(O)=O)C=C1Pathways:Angiogenesis; Immunology/Inflammation; Others; Tyrosine Kinase/AdaptorsReceptor:FGFR1; IL-3; Protein S100-PBoiling pt:Melting pt:>300CSolubility:DMSO: 100 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shishibori T, et al. Biochem J. 1999 Mar 15;338 (Pt 3):583-9.

Note: This product is for in vitro research use only