

Epsin 2 Polyclonal Antibody

Catalog No: #40896

Package Size: #40896-1 50ul #40896-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Epsin 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Epsin 2 Polyclonal Antibody detects endogenous levels of Epsin 2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Epsin 2.
Target Name	Epsin 2
Other Names	EPN2; KIAA1065; Epsin-2; EPS-15-interacting protein 2
Accession No.	Swiss-Prot: O95208NCBI Gene ID: 22905
Uniprot	O95208
GeneID	22905;
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of HuvEc cells using Epsin 2 Polyclonal Antibody

Note: This product is for in vitro research use only