

Pentamidine isethionate

Catalog No: #S1654


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pentamidine isethionate
Brief Description	Inhibitors
Purification	99.83%
Target Name	Others
Calculated MW	592.68
Formulation	C23H36N4O10S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: OCCS(O)(=O)=O.OCCS(O)(=O)=O.NC(=N)C1=CC=C(OCCCCOC2=CC=C(C=C2)C(N)=N)C=C1 Pathways: Others Receptor: Others Boiling pt: 188-194°C Melting pt: Solubility: DMSO: 100mg/ml(168.73 mM); Water: 100mg/ml(168.73 mM); Ethanol:

References

1. Goa KL, Campoli-Richards DM. Drugs. 1987 Mar;33(3):242-58.

Note: This product is for in vitro research use only