

Ambrisentan

Catalog No: #S1655


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ambrisentan
Brief Description	Inhibitors
Purification	99.65%
Target Name	Endothelin Receptor antagonist
Calculated MW	378.42
Formulation	C22H22N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: CC1=CC(=NC(=N1)O[C@H](C(=O)O)C(C2=CC=CC=C2)(C3=CC=CC=C3)OC)C Pathways: GPCR/G

Protein Receptor: ET-ABolig pt: 551.1C? Melting pt: >150C? Solubility: Ethanol: 22 mg/mL (58.13 mM); DMSO: 76 mg/mL (200.83

mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Vatter H, Seifert V. Cardiovasc Drug Rev. 2006 Spring;24(1):63-76.

Note: This product is for in vitro research use only