

Selenomethionine


Catalog No: #S1658

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenomethionine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives
Calculated MW	196.11
Formulation	C ₅ H ₁₁ NO ₂ Se
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: C[Se]CCC(N)C(=O)O
Pathways: Others
Receptor: Amino Acids and Derivatives
Boiling pt: Melting pt: 267-269°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cao J, et al. J Anim Sci Biotechnol. 2014 Oct 2;5(1):46.

Note: This product is for in vitro research use only