

Meglutol


Catalog No: #S1664

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Meglutol
Brief Description	Inhibitors
Purification	98.00%
Target Name	HMG-CoA Reductase inhibitor
Calculated MW	162.14
Formulation	C6H10O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: CC(O)(CC(O)=O)CC(O)=O
Pathways: Metabolism
Receptor: HMG-CoA reductase
Boiling pt: 120°C
Melting pt: 105-108°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Santarelli F, et al. *Ital J Pediatr.* 2013 May 24;39:33.

Note: This product is for in vitro research use only