

Mirabegron

Catalog No: #S1671


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mirabegron
Brief Description	Inhibitors
Purification	99.75%
Target Name	Adrenergic Receptor agonist
Calculated MW	396.51
Formulation	C21H24N4O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: C1=CC=C(C=C1)[C@H](CNCCC2=CC=C(C=C2)NC(=O)CC3=CSC(=N3)N)O Pathways: GPCR/G

Protein Receptor: β 3-adrenoceptor Boiling pt: 138-140°C Melting pt: Solubility: Ethanol: 8 mg/mL (20.17 mM); DMSO: 79 mg/mL (199.23

mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Takasu T, et al. J Pharmacol Exp Ther, 2007, 321(2), 642-647.

Note: This product is for in vitro research use only