


Description

Product Name	Dimesna
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	326.34
Formulation	C ₄ H ₈ Na ₂ O ₆ S ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: [Na+].[Na+].[O-][S](=O)(=O)CCSSCC[S]([O-])(=O)=O
 Pathways: Others
 Receptor: organic anion transporter
 Boiling pt: Melting pt: Solubility: Water: 65 mg/mL (199.17 mM); DMSO: 65 mg/mL (199.17 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Parker AR, et al. Mol, Y Ther, 2010, 9(9), 2558-2567.

Note: This product is for in vitro research use only