

## Acamprosate calcium

Catalog No: #S1673


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Acamprosate calcium
Brief Description	Inhibitors
Purification	99.17%
Target Name	GABA Receptor modulator
Calculated MW	400.48
Formulation	C <sub>10</sub> H <sub>20</sub> CaN <sub>2</sub> O <sub>8</sub> S <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: [Ca++].CC(=O)NCCCS([O-])(=O)=O.CC(=O)NCCCS([O-])(=O)=O  
 Pathways: Membrane transporter/Ion channel; Neuroscience  
 Receptor: GABAR  
 Boling pt: Melting pt: >300C  
 Solubility: Water: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Williams SH. Am Fam Physician. 2005 Nov 1;72(9):1775-80.

Note: This product is for in vitro research use only