

## Nystatin (Fungicidin)

Catalog No: #S1678


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nystatin (Fungicidin)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	926.1
Formulation	C47H75NO17
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CC1OC(OC2CC(O)C(C(O)CC(=O)CC(O)C(O)CCC(O)CC(O)CC(O)CC(=O)OC(C)C(C)C(O)C(C)C=CC=CCCC=CC=C=C=C2)C(O)=O)C(O)C(N)C1O Pathways: Others Receptor: Others Boiling pt: Melting pt: >155°C Solubility: DMSO: 19 mg/mL (20.51 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

- Wallace TL, et al. Antimicrob Agents ChemOther. 1997, 41(10), 2238-2243.

Note: This product is for in vitro research use only