

Tolazoline hydrochloride


Catalog No: #S1680

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tolazoline hydrochloride
Brief Description	Inhibitors
Purification	99.85%
Target Name	Adrenergic Receptor antagonist
Calculated MW	196.68
Formulation	C10H13ClN2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: ClC(C1=NCCN1)C1=CC=CC=C1
Pathways: GPCR/G Protein
Receptor: α -adrenergic receptor
Boiling pt: Melting pt: 172-176°C
Solubility: Ethanol: 39 mg/mL (198.29 mM); Water: 39 mg/mL (198.29 mM); DMSO: 39 mg/mL (198.29 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lei L, et al. J Cardiovasc Pharmacol, 1998, 31(5), 721-733.

Note: This product is for in vitro research use only