

Nitisinone

Catalog No: #S1684


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nitisinone
Brief Description	Inhibitors
Purification	99.54%
Target Name	Others inhibitor
Calculated MW	329.23
Formulation	C ₁₄ H ₁₀ F ₃ NO ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: [O-][N+](=O)C1=C(C=CC(=C1)C(F)(F)F)C(=O)C1C(=O)CCCC1=O Pathways: Others Receptor: 4-hydroxyphenylpyruvate

dioxygenase
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Suwannarat P, et al. Metabolism, 2005. 54(6): p. 719-28.

Note: This product is for in vitro research use only