

## 5Beta-Pregnane-3Alpha,20alpha-Diol

Catalog No: #S1688


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	5Beta-Pregnane-3Alpha,20alpha-Diol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	320.52
Formulation	C21H36O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@@H]([C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC[C@H]4[C@@]3(CC[C@H](C4)O)C)C)O  
 Pathways:Endocrinology/Hormones  
 Receptor:ER  
 Boling pt:Melting pt:Solubility:DMSO: 10 mM  
 Appearance:White Powder  
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Blackwell LF, et al. Hum Reprod. 2013 Dec;28(12):3306-15.

Note: This product is for in vitro research use only