

Enalapril Impurity B

Catalog No: #S1689


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Enalapril Impurity B |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 279.33 |
| Formulation | C ₁₅ H ₂₁ NO ₄ |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
 SMILES: CCOC(=O)[C@H](CCC1=CC=CC=C1)N[C@@H](C)C(=O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Limited solubility
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Chandrakumar NS, et al. J Med Chem. 1992 Jan 24;35(2):223-33.

Note: This product is for in vitro research use only